














Abutment Torque Settings and Screwdriver Type

| Manufacturer | | Torque N-cm | Driver for Titanium Screw |
|---|--|-------------|---------------------------|
|  | ANKLYOS®C?X-C/ 3.5,4.5,5.5,7.0 (non-ind) | 15 | 1.0mm hex |
| | ANKLYOS®C?X-C/ 3.5,4.5,5.5,7.0 (indexed) | 15 | 1.0mm hex |
| | OsseoSpeed™ EV 3.0 - Green | 15 | 0.050" - 1.25mm hex |
| | OsseoSpeed™ EV 3.6 - Purple | 25 | 0.050" - 1.25mm hex |
| | OsseoSpeed™ EV 4.2 - Yellow | 25 | 0.050" - 1.25mm hex |
| | OsseoSpeed™ EV 4.8 - Blue | 25 | 0.050" - 1.25mm hex |
| | OsseoSpeed™ EV 5.4 - Brown | 25 | 0.050" - 1.25mm hex |
| | OsseoSpeed™ Profile EV 4.2 - Yellow | 25 | 0.050" - 1.25mm hex |
| | OsseoSpeed™ Profile EV 4.8 - Blue | 25 | 0.050" - 1.25mm hex |
| | OsseoSpeed™/OsseoSpeed™ TX 3.0 - Yellow | 15 | 0.050" - 1.25mm hex |
| | OsseoSpeed™/OsseoSpeed™ TX 3.5, 4.0 - Aqua | 20 | 0.050" - 1.25mm hex |
| | OsseoSpeed™/OsseoSpeed™ TX 4.5, 5.0 - Lilac | 25 | 0.050" - 1.25mm hex |
|  | NobelActive™ 3.0 | 15 | Unigrip |
| | NobelActive™ NP 3.5 | 35 | Unigrip |
| | NobelActive™ RP 4.3, WP 5.0 | 35 | Unigrip |
| | NobelReplace™ NP 3.5 | 35 | Unigrip |
| | NobelReplace™ NP 4.3 | 35 | Unigrip |
| | NobelReplace™ WP 5.0 | 35 | Unigrip |
| | NobelReplace™ 6.0 | 35 | Unigrip |
| | Branemark System™ NP 3.3 | 20 | 0.048" hex |
| | Branemark System™ RP 3.75, 4.0 | 20 | 0.048" hex |
| Branemark System™ WP 5.0, 5.5 | 20 | 0.048" hex | |
|  | Bone Level 3.3 NC (Narrow CrossFit™) | 35 | SCS |
| | Bone Level 4.1, 4.8 RC (Regular CrossFit™) | 35 | SCS |
| | Standard/Standard Plus 4.8 RN (Regular Neck) | 35 | SCS |
| | Standard/Standard Plus 4.8 WN (Wide Neck) | 35 | SCS |
| | Standard Plus 3.3 NN (Narrow Neck) | 35 | SCS |
|  | Certain™ 3.25, 4.0, 4/3, 5.0, 5/4, 6.0 | 20 | 0.048" hex |
| | Certain™ Prevail™ 3/4/4 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5 | 20 | 0.048" hex |
| | Certain™ XP 4/5, 5/6 | 20 | 0.048" hex |
| | Miniplant™ 3.25 | 20 | 0.048" hex |
| | Osseotite™ 3.75, 4.0, 5.0, 6.0 | 20 | 0.048" hex |
| | Osseotite™ XP 4/5,5/6 | 35 | 0.048" hex |
|  | Internal 3.5, 4.5, 5.7 | 30 | 0.050" - 1.25mm hex |
| | Tapered Internal 3.5, 4.5, 5.7 | 30 | 0.050" - 1.25mm hex |
| | External Implant 3.5 | 30 | 0.050" - 1.25mm hex |
| | External Implant 4.0, 5.0 | 30 | 0.050" - 1.25mm hex |
|  | Screw-line Implant 3.3 | 20 | 0.050" - 1.25mm hex |
| | Screw-Line/Root-line Implant 3.8, 4.3, 5.0, 6.0 | 20 | 0.050" - 1.25mm hex |
|  | Hiossen/Osstem HG mini (3.5) | 30 | 1.2mm hex |
| | Hiossen/Osstem standard 4.0,4.5, 5.0 | 30 | 1.2mm hex |

Abutment Torque Settings and Screwdriver Type

| Manufacturer | | Torque N-cm | Driver for Titanium Screw |
|---|--|-------------|---------------------------|
|  | Imtec® 3M Endure™ | 35 | 0.050" - 1.25mm hex |
| | | | |
|  | Implant Direct Legacy™ | 30 | 0.050" - 1.25mm hex |
| | | | |
|  | MIS® Biocom & Seven | 30 | 0.050" - 1.25mm hex |
| | | | |
|  | OCO Biomedical | 20 | 0.050" - 1.25mm hex |
| | | | |
|  | Thommen SPI® 3.5 | 15 | Four lobe |
| | Thommen SPI® 4.0,4.5,5.0, 6.0 | 20 | Four lobe |
|  | Screw Vent™ 3.3, 3.7, 4.5 | 30 | 0.050" - 1.25mm hex |
| | Screw Vent™ 3.3, 3.7, 4.5 friction fit | 30 | 0.050" - 1.25mm hex |
| | Tapered Screw Vent™ 3.5,4.5, 5.7 | 30 | 0.050" - 1.25mm hex |
| | Tapered Screw Vent™ 3.5, 4.5, 5.7 friction fit | 30 | 0.050" - 1.25mm hex |
| | SwissPlus™ 4.8 | 30 | 0.050" - 1.25mm hex |

CAD/CAM custom titanium and zirconia abutments utilize the same torque and driver settings as the original implant manufacturer