


ATLANTIS™ abutments – Implant compatibility chart



Manufacturer	Implant	ATLANTIS™ Abutment		ATLANTIS™ Crown Abutment		ATLANTIS™ Conus Abutment		Screwdriver
		Titanium Gold-shaded ti	Zirconia	Titanium	Zirconia	Custom	Over-venture	
						Titanium Gold-shaded ti	Titanium	
BioHorizons	INTERNAL							
	Internal 3.5, 4.5, 5.7	x	x	x	x	x	x	0.050" (1.25 mm) Hex screwdriver
	Tapered Internal 3.5, 4.5, 5.7	x	x	x	x	x	x	
	EXTERNAL							
External Implant 3.5	x				x		0.050" (1.25 mm) Hex screwdriver	
External Implant 4.0, 5.0	x	x	x	x	x			
Biomet 3i	INTERNAL							
	Certain™ 3.25, 4.0, 4/3, 5.0, 5/4, 6.0	x	x	x	x	x	x	0.048" Hex screwdriver
	Certain™ Prevail™ 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5	x	x	x	x	x	x	
	Certain™ XP 4/5, 5/6	x	x	x	x	x	x	
	EXTERNAL							
	Miniplant™ 3.25	x	x	x	x	x	x	
Osseofite™ 3.75, 4.0, 5.0, 6.0	x	x	x	x	x	x		
Osseofite™ XP 4/5, 5/6	x	x	x	x	x	x		
Camlog	INTERNAL							
	Screw-line Implant 3.3	x		x		x	x	1.28 mm Hex screwdriver
	Screw-line/Root-line Implant 3.8, 4.3, 5.0, 6.0	x	x	x	x	x	x	
DENTSPLY Implants	INTERNAL							
	ANKYLOS® C/X – C/ 3.5, 4.5, 5.5, 7.0 (non-indexed)	x	x			x		1.0 mm Hex screwdriver
	ANKYLOS® C/X – /X 3.5, 4.5, 5.5, 7.0 (indexed)	x	x			x	x	
	OsseoSpeed™ EV 3.0 – Green	x		x				0.050" Hex screwdriver
	OsseoSpeed™ EV 3.6 – Purple	x	x	x	x	x	x	
	OsseoSpeed™ EV 4.2 – Yellow	x	x	x	x	x	x	
	OsseoSpeed™ EV 4.8 – Blue	x	x	x	x	x	x	
	OsseoSpeed™ EV 5.4 – Brown	x	x	x	x	x	x	
	OsseoSpeed™ Profile EV 4.2 – Yellow	x	x	x		x	x	
	OsseoSpeed™ Profile EV 4.8 – Blue	x	x	x		x	x	
	OsseoSpeed™/OsseoSpeed™ TX 3.0 – Yellow	x		x				
	OsseoSpeed™/OsseoSpeed™ TX 3.5, 4.0 – Aqua	x	x	x	x	x	x	0.050" Hex screwdriver
	OsseoSpeed™/OsseoSpeed™ TX 4.5, 5.0 – Lilac	x	x	x	x	x	x	
	OsseoSpeed™ TX Profile 4.5, 5.0	x		x		x	x	
	OsseoSpeed™ 4.0 – Blue (Previous)	x				x		
XiVE® S 3.0	x		x				1.22 mm Hex screwdriver	
XiVE® S 3.4, 3.8, 4.5, 5.5	x	x	x	x	x	x		
Keystone Dental	INTERNAL							
	PrimaConnex™ SD 3.3, 3.5, Genesis™ 3.8	x	x	x	x	x	x	Keystone Quad Driver
	PrimaConnex™ RD 4.0, 4.1, Genesis™ 4.5	x	x	x	x	x	x	
	PrimaConnex™ WD 5.0, Genesis™ 5.5, 6.5	x	x	x	x	x	x	

Please turn over ➔

Manufacturer	Implant	ATLANTIS™ Abutment		ATLANTIS™ Crown Abutment		ATLANTIS™ Conus Abutment		Screwdriver
		Titanium Gold-shaded ti	Zirconia	Titanium	Zirconia	Custom	Over-enture	
						Titanium Gold-shaded ti	Titanium	
Nobel Biocare	INTERNAL							
	NobelActive™ 3.0	x		x				Nobel Biocare Unigrip
	NobelActive™ NP 3.5	x	x	x	x	x	x	
	NobelActive™ RP 4.3, 5.0	x	x	x	x	x	x	
	NobelReplace™ NP 3.5	x	x	x	x	x	x	
	NobelReplace™ RP 4.3	x	x	x	x	x	x	
	NobelReplace™ WP 5.0	x	x	x	x	x	x	
	NobelReplace™ 6.0	x	x	x	x	x	x	
	NobelReplace™ Conical Connection NP 3.5	x	x	x	x	x	x	
	NobelReplace™ Conical Connection RP 4.3, 5.0	x	x	x	x	x	x	
	EXTERNAL							
	Brånemark System™ NP 3.3	x		x		x	x	Nobel Biocare Unigrip
Brånemark System™ RP 3.75, 4.0	x	x	x	x	x	x		
Brånemark System™ WP 5.0, 5.5	x	x	x	x	x	x		
Straumann	INTERNAL							
	Bone Level 3.3 NC (Narrow CrossFit™)	x	x	x	x	x	x	SCS screwdriver
	Bone Level 4.1, 4.8 RC (Regular CrossFit™)	x	x	x	x	x	x	
	Standard/Standard Plus 4.8 RN (Regular Neck)	x		x		x	x	
	Standard/Standard Plus 4.8 WN (Wide Neck)	x		x		x	x	
	EXTERNAL							
Standard Plus 3.3 NN (Narrow Neck)	x				x		SCS screwdriver	
Zimmer Dental	INTERNAL							
	Screw Vent™ 3.3, 3.7, 4.5	x	x	x	x	x	x	1.25 mm Hex screwdriver
	Screw-Vent™ 3.3, 3.7, 4.5 friction fit	x		x		x	x	
	Tapered Screw-Vent™ 3.5, 4.5, 5.7	x	x	x	x	x	x	
	Tapered Screw-Vent™ 3.5, 4.5, 5.7 friction fit	x		x		x	x	
	SwissPlus™ 4.8	x		x		x	x	

Turnaround times

Turnaround times (in-house business days) for all case sizes		
Titanium	Gold-shaded titanium	Zirconia
2	3	4

Note: Turnaround time listed begins the first business day after receipt of all case information and material required and does not include shipping, which is determined by your choice of delivery method. Turnaround for lab-based scanned cases is one day less than the times shown.

Turnaround times are based on the following:

- Design approvals to process case must be granted by 10 am your local time.
- Lab-scanned cases, including rescans, must be uploaded to ATLANTIS™ WebOrder by 3 pm (EST) in order to save a day of processing compared with the times listed above.
- Material, on-hold and case information issues must be resolved by 8 pm (EST).

Other services

Service	Additional turnaround days
Laboratory services	add 1 day
Provisional restoration	add 5 days
Removable diagnostic preparation	add 2 days

- For the most current implant compatibility, please visit www.dentsplyimplants.com or contact Customer Service at 800-531-3481.
- All orders for ATLANTIS™ abutments are shipped with the corresponding abutment screw needed at no extra charge.

Note: The implant manufacturer's recommended torque should always be used when placing the abutment screw.